


Product Datasheet

Mouse Nppc ELISA Kit (orb565345)

Description	Mouse Nppc ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/ml
Target	Nppc
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q61839
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Nppc(C-
type
natriuretic
peptide)
E...