


Product Datasheet

Mouse Sardh ELISA Kit (orb565446)

Description	Mouse Sardh ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Sardh
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q99LB7
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Sardh(Sarcosine
dehydrogenase,
mit...