


## Product Datasheet

### Mouse Ptpn1 ELISA Kit (orb565467)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse Ptpn1 ELISA Kit |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 0.78-50ng/ml          |
| <b>Target</b>          | Ptpn1                 |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P35821</b>         |
| <b>Sensitivity</b>     | 0.469ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
Ptpn1(Tyrosine-  
protein  
phosphatase...