


Product Datasheet

Mouse Ucn2 ELISA Kit (orb565479)

Description	Mouse Ucn2 ELISA Kit
Reactivity	Mouse
Range	78-5000pg/ml
Target	Ucn2
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q99ML8
Sensitivity	46.9pg/ml
Expiration Date	Please enquire.



Mouse
Ucn2(Urocortin-
2) ELISA Kit