


Product Datasheet

Mouse alpha 2a Adrenergic Receptor ELISA Kit (orb565493)

Description	Mouse alpha 2a Adrenergic Receptor ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Adra2a
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q01338
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Adra2a(Alpha-
2A adrenergic
recepto...