


Product Datasheet

Mouse Rab19 ELISA Kit (orb565499)

Description	Mouse Rab19 ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/ml
Target	Rab19
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P35294
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Mouse
Rab19(Ras-
related
protein
Rab-19) ...