


## Product Datasheet

### Mouse Ephb4 ELISA Kit (orb565506)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse Ephb4 ELISA Kit |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 31.2-2000pg/ml        |
| <b>Target</b>          | Ephb4                 |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P54761</b>         |
| <b>Sensitivity</b>     | 18.75pg/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
Ephb4(Ephrin  
type-B  
receptor 4)  
EL...