


Product Datasheet

Mouse H2-Aa ELISA Kit (orb565509)

Description	Mouse H2-Aa ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/ml
Target	H2-Aa
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P14434
Sensitivity	0.375ng/ml
Expiration Date	Please enquire.



Mouse H2-Aa(H-
2 class II
histocompatibil...