


## Product Datasheet

### Mouse DBP ELISA Kit (orb565569)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse DBP ELISA Kit   |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 3.9-250ng/ml          |
| <b>Target</b>          | DBP/GC                |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P21614</b>         |
| <b>Sensitivity</b>     | 2.34ng/ml             |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
DBP/GC(Vitamin  
D-binding  
protein) ...