


## Product Datasheet

### Mouse Trim72 ELISA Kit (orb565630)

|                        |                        |
|------------------------|------------------------|
| <b>Description</b>     | Mouse Trim72 ELISA Kit |
| <b>Reactivity</b>      | Mouse                  |
| <b>Range</b>           | 78-5000pg/ml           |
| <b>Target</b>          | Trim72                 |
| <b>Storage</b>         | 2-8 °C for 6 months    |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>Q1XH17</b>          |
| <b>Sensitivity</b>     | 46.9pg/ml              |
| <b>Expiration Date</b> | Please enquire.        |



Mouse  
Trim72(Tripartite  
motif-  
containing...