


Product Datasheet

Mouse Prl8a8 ELISA Kit (orb565658)

Description	Mouse Prl8a8 ELISA Kit
Reactivity	Mouse
Range	78-5000pg/ml
Target	Prl8a8
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q9DAS4
Sensitivity	46.9pg/ml
Expiration Date	Please enquire.



Mouse
Prl8a8(Prolactin-
8A8) ELISA Kit