


Product Datasheet

Mouse Slco1a1 ELISA Kit (orb565677)

Description	Mouse Slco1a1 ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/ml
Target	Slco1a1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q9QXZ6
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Mouse
Slco1a1(Solute
carrier organic
ani...