


Product Datasheet

Mouse ANP ELISA Kit (orb565693)

Description	Mouse ANP ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	ANP
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P05125
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
ANP(Atrial
Natriuretic
Peptide)
EL...