


## Product Datasheet

### Mouse COL1 ELISA Kit (orb565778)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse COL1 ELISA Kit  |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 0.313-20ng/ml         |
| <b>Target</b>          | COL1                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P11087</b>         |
| <b>Sensitivity</b>     | 0.188ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
COL1(Collagen  
Type I) ELISA  
Kit