


## Product Datasheet

### Mouse CX26 ELISA Kit (orb565799)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse CX26 ELISA Kit  |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 19.531-1250pg/ml      |
| <b>Target</b>          | CX26                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q00977</b>         |
| <b>Sensitivity</b>     | 11.719pg/ml           |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
CX26(Connexin  
26) ELISA Kit