


Product Datasheet

Human Asprosin protein (orb5658)

Description	Recombinant of human Asprosin protein (N-8His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Protein Sequence	Ser2732-His2871
Source	Human Cells
Uniprot ID	P35555
Expiration Date	6 months from date of receipt.



Recombinant
Human
Fibrillin-
1/Asprosin
(...)