


## Product Datasheet

### Mouse CX3CR1 ELISA Kit (orb565800)

|                        |                        |
|------------------------|------------------------|
| <b>Description</b>     | Mouse CX3CR1 ELISA Kit |
| <b>Reactivity</b>      | Mouse                  |
| <b>Range</b>           | 0.156-10ng/ml          |
| <b>Target</b>          | CX3CR1                 |
| <b>Storage</b>         | 2-8 °C for 6 months    |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>Q9Z0D9</b>          |
| <b>Sensitivity</b>     | 0.094ng/ml             |
| <b>Expiration Date</b> | Please enquire.        |



Mouse  
CX3CR1(Chemokine  
C-X3-C-Motif Rece...