


Product Datasheet

Mouse DPD ELISA Kit (orb565822)

Description	Mouse DPD ELISA Kit
Reactivity	Mouse
Range	3.125-200ng/ml
Target	DPD
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q8CHR6
Sensitivity	1.875ng/ml
Expiration Date	Please enquire.



Mouse
DPD(Deoxyipyridinoline)
ELISA Kit