


## Product Datasheet

### Mouse ERN1 ELISA Kit (orb565834)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse ERN1 ELISA Kit  |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 39.063-2500pg/ml      |
| <b>Target</b>          | ERN1                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q9EQY0</b>         |
| <b>Sensitivity</b>     | 23.438pg/ml           |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
ERN1(Endoplasmic  
Reticulum To  
Nucl...