


## Product Datasheet

### Mouse FGB ELISA Kit (orb565846)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse FGB ELISA Kit   |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 3.906-250ng/ml        |
| <b>Target</b>          | FGB                   |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q8K0E8</b>         |
| <b>Sensitivity</b>     | 2.344ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
FGB(Fibrinogen  
Beta chain)  
ELISA K...