


Product Datasheet

Mouse GSTalpha1 ELISA Kit (orb565907)

Description	Mouse GSTalpha1(Glutathione S Transferase Alpha 1) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	GSTα1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P13745
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
GSTα1(Glutathione
S Transferase A...