


## Product Datasheet

### Mouse Gzms-A ELISA Kit (orb565911)

|                        |                        |
|------------------------|------------------------|
| <b>Description</b>     | Mouse Gzms-A ELISA Kit |
| <b>Reactivity</b>      | Mouse                  |
| <b>Range</b>           | 0.156-10ng/ml          |
| <b>Target</b>          | Gzms-A                 |
| <b>Storage</b>         | 2-8 °C for 6 months    |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>P11032</b>          |
| <b>Sensitivity</b>     | 0.094ng/ml             |
| <b>Expiration Date</b> | Please enquire.        |



Mouse  
Gzms-  
A(granzyme  
A) ELISA Kit