


## Product Datasheet

### Mouse HDL ELISA Kit (orb565920)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse HDL ELISA Kit   |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 0.313-20nmol/ml       |
| <b>Target</b>          | HDL                   |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Sensitivity</b>     | 0.188nmol/ml          |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
HDL(High  
Density  
Lipoprotein)  
ELIS...