


## Product Datasheet

### Mouse NOX4 ELISA Kit (orb566020)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse NOX4 ELISA Kit  |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 0.156-10ng/ml         |
| <b>Target</b>          | NOX4                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q9JH18</b>         |
| <b>Sensitivity</b>     | 0.094ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
NOX4(Nicotinamide  
Adenine Dinucleo...