


Product Datasheet

Mouse PDGFsR-alpha ELISA Kit (orb566050)

Description	Mouse PDGFsR-alpha(Platelet-Derived Growth Factor Soluble Receptor alpha) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	PDGFsR- α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P26618
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
PDGFsR- α (Platelet-Derived Growth ...