


Product Datasheet

Mouse PPAR-gamma ELISA Kit (orb566071)

Description	Mouse PPAR-gamma(Peroxisome Proliferator-activated receptor gamma) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	PPAR-γ
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P37238
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse PPAR-
γ(Peroxisome
Proliferator-
ac...