


Product Datasheet

Mouse SOD3 ELISA Kit (orb566124)

Description	Mouse SOD3 ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	SOD3
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	O09164
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
SOD3(Superoxide
Dismutase 3,
Extra...