


## Product Datasheet

### Mouse TFPI2 ELISA Kit (orb566147)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse TFPI2 ELISA Kit |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 0.156-10ng/ml         |
| <b>Target</b>          | TFPI2                 |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>O35536</b>         |
| <b>Sensitivity</b>     | 0.094ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
TFPI2(Tissue  
Factor  
Pathway  
Inhibi...