


Product Datasheet

Mouse alpha1-AGP ELISA Kit (orb566183)

Description	Mouse alpha1-AGP(alpha1-Acid glycoprotein) ELISA Kit
Reactivity	Mouse
Range	6.25-400ng/ml
Target	α 1-AGP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q60590
Sensitivity	3.75ng/ml
Expiration Date	Please enquire.



Mouse α 1-AGP(α 1-Acid glycoprotein) ELISA...