



Product Datasheet

Acetyl-Histone H4 (Lys8) antibody (orb5662)

Description	Acetyl-Histone H4 (Lys8) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated Synthesised acetylpeptide derived from human Histone H4 around the acetylation site of Lys8 (GG(Ac-K)GL)
Target	Histone H4 (Acetyl K8)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11 kDa
Uniprot ID	P62805
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1ug/test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human cell; Lane...



ICC analysis of HeLa cell using Histone ...