


Product Datasheet

Mouse Heparan Sulfate ELISA Kit (orb566251)

Description	Mouse Heparan Sulfate ELISA Kit
Reactivity	Mouse
Range	31.25-2000ng/ml
Target	Heparan Sulfate
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sensitivity	18.75ng/ml
Expiration Date	Please enquire.



Mouse
HS(Heparan
Sulfate)
ELISA Kit