


Product Datasheet

Mouse Prl3b1 ELISA Kit (orb566253)

Description	Mouse Prl3b1 ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Prl3b1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P09586
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Prl3b1(Prolactin
3B1) ELISA Kit