


Product Datasheet

Mouse Ang1-7 ELISA Kit (orb566285)

Description	Mouse Ang1-7 ELISA Kit
Reactivity	Mouse
Range	7.8-500pg/ml
Target	Ang1-7
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P11859
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Mouse Ang1-7(Angiotensin 1-7) ELISA Kit