


Product Datasheet

Mouse CAMP ELISA Kit (orb566288)

Description	Mouse CAMP ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	CAMP
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Sandwich ELISA, Double Antibody
Uniprot ID	P51437
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
CAMP(Cathelicidin
Antimicrobial Pe...