


## Product Datasheet

### Mouse AFP-L3 ELISA Kit (orb566316)

|                        |                        |
|------------------------|------------------------|
| <b>Description</b>     | Mouse AFP-L3 ELISA Kit |
| <b>Reactivity</b>      | Mouse                  |
| <b>Range</b>           | 0.313-20ng/ml          |
| <b>Target</b>          | AFP-L3                 |
| <b>Storage</b>         | 2-8 °C for 6 months    |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>P02772</b>          |
| <b>Sensitivity</b>     | 0.188ng/ml             |
| <b>Expiration Date</b> | Please enquire.        |



Mouse  
AFP-  
L3(Alpha-  
Fetoprotein  
L3) ELISA...