


Product Datasheet

Mouse alpha 1 Antitrypsin ELISA Kit (orb566320)

Description	Mouse alpha 1 Antitrypsin ELISA Kit
Reactivity	Mouse
Range	31.25-2000ng/ml
Target	AAT
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P07758
Sensitivity	18.75ng/ml
Expiration Date	Please enquire.



Mouse AAT
(Alpha 1
Antitrypsin)
ELISA Ki...