


## Product Datasheet

### Mouse Rbks ELISA Kit (orb566352)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Mouse Rbks ELISA Kit  |
| <b>Reactivity</b>      | Mouse                 |
| <b>Range</b>           | 78.125-5000pg/ml      |
| <b>Target</b>          | Rbks                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q8R1Q9</b>         |
| <b>Sensitivity</b>     | 46.875pg/ml           |
| <b>Expiration Date</b> | Please enquire.       |



Mouse  
Rbks(Ribokinase)  
ELISA Kit