


Product Datasheet

Monkey TNFRSF1A ELISA Kit (orb566421)

Description	Monkey TNFRSF1A ELISA Kit
Reactivity	Monkey
Range	15.625-1000pg/ml
Target	TNFRSF1A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	A0A2K6U4D0
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Monkey
TNFRSF1A(Tumor
Necrosis Factor
Re...