


Product Datasheet

Monkey TNF-alpha ELISA Kit (orb566422)

Description	Monkey TNF-alpha(Tumor Necrosis Factor Alpha) ELISA Kit
Reactivity	Monkey
Range	31.25-2000pg/ml
Target	TNF- α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P79337
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Monkey
TNF-
 α (Tumor
Necrosis
Factor
Alpha...)