


Product Datasheet

Monkey alphaFP ELISA Kit (orb566425)

Description	Monkey alphaFP(Alpha-Fetoprotein) ELISA Kit
Reactivity	Monkey
Range	1.563-100ng/ml
Target	α FP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	A0A2K5U7T6
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Monkey
 α FP(Alpha-
Fetoprotein)
ELISA Kit