


Product Datasheet

Monkey GNRH ELISA Kit (orb566453)

Description	Monkey GNRH (Gonadotropin-releasing hormone) ELISA Kit
Reactivity	Monkey
Range	15.625-1000nmol/ml
Target	GNRH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	A0A2K5VL54
Sensitivity	9.375nmol/ml
Expiration Date	Please enquire.



Monkey
GNRH(Gonadotropin-releasing hormo...