


## Product Datasheet

### Porcine ET-1 ELISA Kit (orb566492)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Porcine ET-1 ELISA Kit                 |
| <b>Reactivity</b>        | Porcine                                |
| <b>Range</b>             | 1.25-80pg/ml                           |
| <b>Target</b>            | ET-1                                   |
| <b>Storage</b>           | 2-8 °C for 6 months                    |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Competitive ELISA, Coated with Antigen |
| <b>Uniprot ID</b>        | <b>P09558</b>                          |
| <b>Sensitivity</b>       | 0.625pg/ml                             |
| <b>Expiration Date</b>   | Please enquire.                        |



Porcine ET-1(Endothelin 1) ELISA Kit