


## Product Datasheet

### Porcine LH ELISA Kit (orb566541)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Porcine LH ELISA Kit                   |
| <b>Reactivity</b>        | Porcine                                |
| <b>Range</b>             | 0.781-50ng/ml                          |
| <b>Target</b>            | LH                                     |
| <b>Storage</b>           | 2-8 °C for 6 months                    |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Competitive ELISA, Coated with Antigen |
| <b>Uniprot ID</b>        | <b>P01232</b>                          |
| <b>Sensitivity</b>       | 0.469ng/ml                             |
| <b>Expiration Date</b>   | Please enquire.                        |



Porcine  
LH(Luteinizing  
Hormone)  
ELISA Ki...