


Product Datasheet

Porcine VEGFR1 ELISA Kit (orb566587)

Description	Porcine VEGFR1/Flt1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/ml
Target	VEGFR1/Flt1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	F1RSU5
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Porcine
VEGFR1/Flt1(Vascuoar
Endothelial...