


Product Datasheet

Porcine alphaFP ELISA Kit (orb566590)

Description	Porcine alphaFP(Alpha-Fetoprotein) ELISA Kit
Reactivity	Porcine
Range	0.313-20ng/ml
Target	α FP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8MJ76
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Porcine
 α FP(Alpha-
Fetoprotein)
ELISA Ki...