


Product Datasheet

Porcine alpha1-AGP ELISA Kit (orb566624)

Description	Porcine alpha1-AGP(alpha1-Acid Glycoprotein) ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/ml
Target	α 1-AGP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	F1SN68
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Porcine α 1-AGP(α 1-Acid Glycoprotein) E...