


Product Datasheet

Rat Endothelin 1 ELISA Kit (orb566665)

Description	Rat Endothelin 1 ELISA Kit
Reactivity	Rat
Range	1.25-80pg/ml
Target	ET-1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P22388
Sensitivity	0.625pg/ml
Expiration Date	Please enquire.



Rat ET-1(Endothelin 1) ELISA Kit