


## Product Datasheet

### Rat IGFBP-3 ELISA Kit (orb566714)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Rat IGFBP-3 ELISA Kit |
| <b>Reactivity</b>      | Rat                   |
| <b>Range</b>           | 78.125-5000pg/ml      |
| <b>Target</b>          | IGFBP-3               |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P15473</b>         |
| <b>Sensitivity</b>     | 46.875pg/ml           |
| <b>Expiration Date</b> | Please enquire.       |



Rat  
IGFBP-  
3(Insulin-  
like  
growth  
factor-b...