


Product Datasheet

Rat HIF-1alpha ELISA Kit (orb566761)

Description	Rat HIF-1alpha(Hypoxia-inducible factor 1-alpha) ELISA Kit
Reactivity	Rat
Range	7.813-500pg/ml
Target	HIF-1 α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O35800
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Rat HIF-1 α (Hypoxia-inducible factor 1-a...