


## Product Datasheet

### Rat Npr2 ELISA Kit (orb566953)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Rat Npr2 ELISA Kit    |
| <b>Reactivity</b>      | Rat                   |
| <b>Range</b>           | 78-5000pg/ml          |
| <b>Target</b>          | Npr2                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P16067</b>         |
| <b>Sensitivity</b>     | 46.9pg/ml             |
| <b>Expiration Date</b> | Please enquire.       |



Rat  
Npr2(Atrial  
natriuretic  
peptide  
rece...